

BIO-PLUS LIGHT
REGENflex

**CROSS-LINKED
HYALURONIC ACID**
40 mg / 1,6 ml

REGENYAL
LABORATORIES

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DESCRIPTION OF FUNCTIONAL PROPERTIES

REGENFLEX® BIO-PLUS LIGHT is a buffered physiological solution, containing 2,5% of cross-linked hyaluronic acid with significant viscous and elastic properties.

The sodium hyaluronate is a fundamental component of the synovial fluid, being responsible for the particular viscous and elastic properties of the latter.

REGENFLEX® BIO-PLUS LIGHT

1 pre-filled mono-use syringe (40 mg of reticulated hyaluronic acid in 1,6 ml buffered physiological solution).

Moist heat sterilized syringe.

COMPOSITION

Main ingredient: cross-linked hyaluronic acid.

OTHER INGREDIENTS:

Sodium chloride, phosphate buffer, water for injectable solutions.

INDICATIONS

REGENFLEX® BIO-PLUS LIGHT is a substitute for the synovial fluid that can restore physiological and rheological properties of compromised joints, when there are pain or reduced mobility caused by degenerative or post-traumatic diseases.

DOSAGE AND INSTRUCTIONS FOR USE

REGENFLEX® BIO-PLUS LIGHT: Inject 1,6 ml of the product into the joint every 3 months, unless otherwise decided by the physician; each treatment should be tailored to the individual patient.

Depending on the physician's opinion and the patient's health, it is possible to contemporarily treat several joints.

The volume scale on the syringes is merely indicative.

MODE OF USE

Before injecting REGENFLEX® BIO-PLUS LIGHT, remove by aspiration any eventual articular effusion. Carefully remove the cap from the syringe in order to avoid any contact with the opening. Insert an appropriate needle (from 18 to 22G) twisting it on to ensure that it is air tight and to prevent leakage of the solution.

Inject the medical device only into the synovial space, at room temperature and in complete asepsis.

WARNINGS AND CONTRAINDICATIONS

Do not use after the expiry date printed on the box.

Do not use if the package is open or damaged.

Injections must be performed on healthy skin.

Do not inject into blood vessels.

Do not inject into the knee-joint when there are venous or lymphatic stasis of the limb.

Do not inject outside the joint cavity, into the synovial tissue or into the joint capsule.

Do not administer in case of excessive intra-articular effusions.

REGENFLEX® BIO-PLUS LIGHT must not be administered in case of infections, inflammations or skin diseases in the concerned area.

After the intra-articular injection the patient should be advised to avoid any intense physical activity and to take up normal activities only after a few days.

The content of the pre-filled syringe is sterile.

The external surface of the syringe is not sterile.

Do not sterilize anew. The device is mono-use.

Store at temperatures lower than 28 °C and away from heat.

Do not freeze.

Once opened the device must be used immediately and disposed after use.

Keep out of the reach of children.

SIDE EFFECTS

Extra-articular injection may cause local side effects.

The injection area might become painful, hot, reddened or swollen. When this is the case, it can be useful to apply ice to the treated joint. Such symptoms are normal reactions and disappear within a short time.

The physician must ensure that the patient inform them of any eventual side effect that may occur after the treatment.

In rare cases, the product can cause allergic reactions.

REGENFLEX® BIO-PLUS LIGHT IS AN ON PRESCRIPTION DEVICE.

INTRA-ARTICULAR INJECTION MUST ONLY BE PERFORMED BY A PHYSICIAN.

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REGENYAL LABORATORIES Srl